



Accuracy Matters When Quantitative, Manually-pipetted Assays Graduate to Automation – A Story in Diagnosing and Troubleshooting

Keith J. Albert (ARTEL*) and Lisa M. Knapp (Beckman Coulter, 7451 Winton Drive, Indianapolis, IN 46268)

Abstract

It is often the case that assays are initially performed on the benchtop using handheld pipettes before they graduate, or transfer, to an automated liquid handler. Automating a manual method may take time and patience, but automation will help lower costs, increase throughput, and potentially avoid errors associated with a manual method. During the transfer process, however, the manual assay should be directly compared to the automated assay for consistencies in pipetting performance. What if the volumes of reagents pipetted to and from the assays are only compared for precision and not accuracy? An undetected variability in accuracy will impact the integrity of the assay as the automation process continues. Even if the manual method and automated method both have highly precise liquid delivery, they might not be pipetting the same volume of reagent to the assay. Accuracy information is critical to determine the deviation of the dispensed volumes from the target volume when transferring or scaling up an assay. Validating the liquid delivery steps for each assay using a robust measurement method should be implemented to uncover discrepancies in pipetting performance. This presentation discusses the importance of knowing both accuracy and precision information when a manual method is transferred to a Biomek® liquid handler. By using a standardized volume verification methodology, it was determined that the rate-limiting reagent was not being accurately pipetted between the manual and automated methods per the protocol, and the automation was not to blame.

Introduction

The need to ensure quality in a laboratory process has become increasingly important, especially as it relates to the accurate volume transfers within microplate-based assays (ref 1). Measuring and knowing the exact volume transferred will inherently lead to accurate and precise analysis of the experiment, i.e., the results can be trusted. Many assays, such as the highlighted real time polymerase chain reaction (RT-PCR) assay, depend on accurate volume delivery with either a handheld pipette (benchtop assay) or a liquid handler (automated assay). Because the concentrations of components in assays are volume dependent, inaccurate volume transfers will directly impact assay results. In our story, the RT-PCR assay was being transferred from a manual pipetted method to an automated liquid handler (Biomek FX). The technician was observing acceptable results with the manual assay, but as the method graduated, or transferred, to the Biomek, the results could not be repeated and were deemed unacceptable, i.e., the liquid handler was producing errant results for the same RT-PCR assay that worked with a handheld pipette. The liquid handler was blamed for performing poorly and a service call was initiated. A senior

Introduction (cont'd)

applications scientist was called-in to diagnose, troubleshoot and correct errors in the automated liquid handling method. The field service engineer employed the MVS® Multichannel Verification System (Figure 1-left), which is a standardized volume measurement platform, to quantify the volume transfer performance of the Biomek FX (Figure 1-right).



Figure 1. The ARTEL MVS® Multichannel Verification System (left) was used to assess the volume transfer performance for the Beckman Coulter Biomek FX (right) for the RT-PCR assay.

Accuracy required for the RT-PCR assay

This RT-PCR assay was primer-limited, so the volume transfer of 15 µL master mix to the assay was a critical step and therefore needed to be accurate.

RT – PCR assay

- Assay Plate
 - 10 µL Sample – In Excess; not a critical volume transfer step
 - 15 µL Master Mix – **CRITICAL STEP; accuracy required**
- Place in thermocycler for 1.5 hours
- Results obtained

The applications scientist used the MVS to characterize the Biomek's performance when transferring the 15-µL target volume. As discussed in reference 2 and shown in the schematic in Figure 2, the MVS dye-based sample solutions are dispensed into the microplate by the liquid handler, and the target volume is quantified well-by-well and tip-by-tip using dual-eye radiometric photometry with a foundation based on the Beer-Lambert Law. Immediately after the absorption measurements are collected, the MVS is used to determine the volume dispensed in each well of a microplate (ref 2).




Figure 2. The volume transfer performance for automated liquid handling and hand pipetting can be assessed with the radiometric photometry implemented by the MVS. A target volume of sample solution is dispensed into a microplate well before, or after, the non-quantitative addition of diluent (ref 2). The solutions are mixed and the absorbance values are used to determine the well-by-well volume values using a series of equations. The MVS is a standardized, universal platform for measuring both accuracy and precision of volume delivery in 96-well or 384-well plates for volumes between 10 nL – 200 µL. Such measurements allow for directly comparing pipetting equipment regardless of tip number, model, make, manufacturer or location.

“As Found” Performance Assessment

The Biomek's “as found” volume transfer performance was determined with the MVS for the 15-µL target volume of aqueous solution in a 96-well plate. The initial Scaling Factor and Offset for the Biomek were set at 1 and 0, respectively.

This partial MVS Output Report for “as found” data shows that the Biomek's accuracy is out of tolerance with an average transferred volume of 14.03 µL.

Overall test statistics flagged for failed accuracy

Tip	Volume (µL)	Mean	Std. Dev.	CV (%)	Rel. Inacc. (%)	CV Flag	Rel. Inacc. Flag
1	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
2	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
3	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
4	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
5	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
6	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
7	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
8	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
9	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
10	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
11	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
12	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
13	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
14	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
15	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
16	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
17	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
18	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
19	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
20	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
21	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
22	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
23	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
24	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
25	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
26	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
27	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
28	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
29	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
30	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
31	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
32	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
33	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
34	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
35	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
36	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
37	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
38	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
39	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
40	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
41	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
42	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
43	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
44	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
45	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
46	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
47	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
48	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
49	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
50	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
51	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
52	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
53	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
54	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
55	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
56	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
57	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
58	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
59	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
60	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
61	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
62	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
63	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
64	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
65	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
66	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
67	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
68	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
69	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
70	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
71	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
72	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
73	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
74	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
75	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
76	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
77	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
78	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
79	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
80	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
81	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
82	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
83	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
84	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
85	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
86	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
87	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
88	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
89	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
90	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
91	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
92	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
93	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
94	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
95	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail
96	15.00	14.03	0.07	0.51%	-4.47%	Pass	Fail

Figure 3. A partial MVS output report for the initial, “as found”, performance assessment shows highlighted (flagged) tips and wells that perform outside of the assay specifications. Assay specifications for relative inaccuracy and CV were set at 3% and 5%, respectively. Though the CV for all tips was within specification (at 1.14%), the relative inaccuracy for all tips was falling at -6.47% with an average dispensed volume for the 8 tips and 3 reps of 14.03 µL.

Adjusting the Software for Accuracy

In the Biomek software, the Scaling Factor and Offset were adjusted to 1.045 and 0.354, respectively and the 15-µL target volume was re-tested to show that the instrument was dispensing within the assay's specifications (Figure 4).

After adjustments, the Biomek's accuracy and precision values are within assay specifications with a measured 15.01-µL average volume dispensed.

Tip	Volume (µL)	Mean	Std. Dev.	CV (%)	Rel. Inacc. (%)	CV Flag	Rel. Inacc. Flag
1	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
2	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
3	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
4	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
5	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
6	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
7	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
8	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
9	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
10	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
11	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
12	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
13	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
14	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
15	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
16	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
17	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
18	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
19	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
20	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
21	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
22	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
23	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
24	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
25	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
26	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
27	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
28	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
29	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
30	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
31	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
32	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
33	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
34	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
35	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
36	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
37	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
38	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
39	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
40	15.00	15.01	0.07	0.47%	0.67%	Pass	Pass
41	15.00	15.01	0.07	0.47%	0.67%		