



**ARTELM**

***SLAS2019 Tutorial: Coupling Assay Design And  
Process Optimization Toward Minimizing Variability***

***Nathaniel Hentz, PhD  
05 February; 12:30-1:15P***


***Acknowledgment: Becky Kitchener, PhD***

# Tutorial Agenda

- Assay variability
- Liquid handler optimization: current practices
- Process optimization – What? Why? How?
- Case study: model assay & findings
- Conclusions

# OVERVIEW: ASSAY VARIABILITY

# Parameters that Effect Assay Variability



# Assay Validation in a High-Throughput World



# So, What Does A “Good” Assay Look Like?



Let’s define a couple of parameters:

- First, let’s consider the assay variability
- Second, let’s consider the assay window


$$Z = 1 - \left( \frac{3\sigma_{\max} + 3\sigma_{\min}}{|\mu_{\max} - \mu_{\min}|} \right)$$

Z-Factor Value	Assay Quality
1	Ideal assay
$1 > Z \geq 0.5$	Excellent assay
$0.5 > Z > 0$	Marginal assay
0	Yes/no assay
$< 0$	Assay not useful

# The Result



# The Result





# OVERVIEW: LIQUID HANDLER OPTIMIZATION

# Liquid Dispense Technologies

Dispense Technology	Attributes
<b>Air displacement</b>	<ul style="list-style-type: none"><li>• Problematic for volatile liquids</li><li>• Possible cross contamination</li><li>• Wide range of volumes</li></ul>
<b>Positive displacement</b>	<ul style="list-style-type: none"><li>• Useful for volatile solvents</li><li>• Cross contamination possible</li><li>• Wide range of volumes</li></ul>
<b>Droplet (acoustic)</b>	<ul style="list-style-type: none"><li>• Non-contact</li><li>• Useful for small volumes (pL – nL)</li></ul>
<b>Droplet (solenoid/inkjet)</b>	<ul style="list-style-type: none"><li>• Useful for small volumes (nL)</li><li>• Sensitive to fluid types</li></ul>
<b>Capillary (pintool)</b>	<ul style="list-style-type: none"><li>• Useful for small volumes</li><li>• Direct contact with sample (contamination)</li><li>• Sensitive to fluid type</li></ul>
<b>Peristaltic</b>	<ul style="list-style-type: none"><li>• Useful for bulk dispense;</li><li>• More frequent calibrations needed</li></ul>

## Parameters that Effect Volume Transfer


- **Hardware/labware** – plate format, tip types, tubing type
- **ALH settings** – target (or off-set) volume, single/multi dispense, aspirate/dispense rate, aspirate/dispense height, liquid class, pre and post air gaps, accuracy/precision of volume transfer, transfer speed/time delays, on-board mixing
- **Assay** – reagent mixing, incubation, centrifugation, number of liquid transfers, wash steps, dilution steps, serial dilutions
- **Sample** – viscosity, density, surface tension, temperature, volatility
- **Environment** – temperature, humidity, light, vibration

# Typical Liquid Handler Optimization

- Usually performed as a stand-alone activity
  - Precision is always checked, but accuracy is not as easy
- Several methods for volume verification
  - In-house (e.g., fluorescence, gravimetric, absorbance, etc.)
  - Commercial (e.g., dual-dye spectrophotometry)
- Volume verification is typically performed with ideal solutions
- Liquid handler is certified, calibrated, or repaired (if necessary)
- Then....someone programs ALH for assay use
  - Default method
  - Specific to basic assay requirements

# Artel MVS: A Useful Tool

- Employs a dual-dye, dual-wavelength, ratiometric absorbance-based measurement method for calculating the dispense volume.
- How it works: dyes of known concentration are dispensed into well-characterized microtiter plate. The plate is mixed on a plate shaker to ensure solution homogeneity. Absorbance readings are taken at 520 nm and 730 nm.



# ADOPTING A PROCESS OPTIMIZATION APPROACH

# Sources of Assay Variability

*What can we control or optimize?*

## Biology

- Diffusion
- Binding equilibrium
- Steric hindrance
- Cell population diversity
- Protein activity

## Environment

- ✓ Temperature
- ✓ Light
- ✓ Humidity
- People

## Instrumentation

- ✓ Liquid handlers
- ✓ Detectors
- ✓ Mixers
- ✓ Incubators
- ✓ Centrifuges

## Consumables

- ✓ Reagents
- ✓ Tips
- ✓ Plate type
- ✓ Plate seals

# Artel MVS is More Than a Calibration Tool





# CASE STUDY

# Assessing Effects of Liquid Handling on the Assay


- **Phase I:** Using a well-characterized model assay, develop and optimize the assay platform at bench scale.
- **Phase II:** Perform method transfer to ALH: examine effects of automated liquid handling parameters on the same assay.
- **Phase III:** “Deconstruct” the assay: decide which parameters, when altered, significantly vary assay outcome.

# Streptavidin: Biotin-FI Assay Principle

Streptavidin (SA) is a tetravalent biotin-binding protein that is isolated from *Streptomyces avidinii* and has a mass of 60.0 kDa. SA has a very high affinity for biotin ( $K_d = 10^{-14}$  to  $10^{-15}$  M).

## Assay Set-up

- PBS, PBS+0.1%BSA, and PBS+0.1% glycerol
- Add 25  $\mu$ L of biotin-4-fluorescein (FI) to black 96w plate
- Add 25  $\mu$ L of inhibitor
- Add 25  $\mu$ L of streptavidin (SA)
- Incubate for 30 min at room temperature
- Read fluorescence: Ex = 485 nm and Em = 515 nm




1. Waner, MJ; Mascotti DP. *Journal of Biochemical and Biophysical Methods* 70(6), 2008, 873-877.
2. Ebner, A; Marek, M; Kaiser, K; Kada, G; Hahn, CD; Lackner, B; Gruber, HJ. *Methods in Molecular Biology*, 418, 2008, 73-88.

# Assay Development Parameters Evaluated

Parameter	Manual
Solution stability	X
Plate type	X
Background buffer	X
Inhibitor	X
Inhibitor conc.	X
Incubation time	X
Mixing	X
Pipette calibration	
Asp/Disp rates	
Fluid exit rate	
Air gaps	

# Effect of Mixing on Assay Variability





*Mixing was conducted by aspirating and dispensing 3 cycles after the third reagent was added. The plates were then incubated for 30 minutes at room temperature.*

## Key Takeaway

Buffer has an impact on mixing; Mixing improved BSA and glycerol buffers

PBS = phosphate buffered saline / BSA = PBS + 0.1% bovine serum albumin / GL = PBS + 0.1% glycerol

# Effect of Buffer Type on Assay Performance




## Key Takeaway

Manual mixing can behave differently than automated mixing, especially depending on assay buffer ingredients.

# Effect of Source Plates on Automated Liquid Handler Verification

**A**



**B**


Artel MVS plates containing reagent from non-binding (A) and untreated (B) 96w black source plates.

MVS precision: 6.3% for non-binding and 1.2% for untreated

## Key Takeaway

Source plate affected the destination plate – observed by visual inspection

# Effect of Plate Type on Assay Variability




- Black, 96-well plates used:
- Corning #3650 (Non-binding)
  - Corning #3915 (Untreated)

**Key Takeaway**  
When subjected to the assay, the uncoated plate yielded slightly lower variability.

PBS = phosphate buffered saline / BSA = PBS + 0.1% bovine serum albumin / GL = PBS + 0.1% glycerol




# Effect of Reagent Temperature on Assay Variability




“Cold reagents” were stored at 4°C until use. The reagents were not maintained at 4°C during pipetting.

# Effect of Aspirate (ASP) Speed on Assay Performance




PBS = phosphate buffered saline / BSA = PBS + 0.1% bovine serum albumin / GL = PBS + 0.1% glycerol

# Effect of Dispense (DISP) Speed on Assay Performance



PBS = phosphate buffered saline / BSA = PBS + 0.1% bovine serum albumin / GL = PBS + 0.1% glycerol


# Effect of Air Gap on Inhibitor Potency



**Key Takeaway**  
 Air gap affects the assay containing PBS and BSA buffers with respect to potency and variability

PBS = phosphate buffered saline / BSA = PBS + 0.1% bovine serum albumin / GL = PBS + 0.1% glycerol

# Effect of Air Gap on Assay Variability



PBS = phosphate buffered saline / BSA = PBS + 0.1% bovine serum albumin / GL = PBS + 0.1% glycerol

## Study Summary


- Four assay parameters were identified which required optimization both during development at bench-scale *and* on the ALH.
  - Assay buffer selection
  - Mixing
  - Reagent stability
  - Plate type
- Certain ALH parameters were dependent on buffer type. Not all ALH parameters will effect an assay.

# Conclusions

- Performing assay optimization **AND** liquid handler optimization **together as a whole process** reduces potential for error introduction and prolonged/difficult method transfer.
- Evaluate critical liquid handling parameters and potential sources of variability at bench scale and on the ALH for **each new assay**.

**Assay optimization or LH qualification alone isn't enough.**

# Putting It Together...



**Assay results are impacted by liquid handling variability**





**ARTEL**

*Trust Your Results*